

RNA isolation from Tempus blood RNA tubes by Versagene RNA purification kit

The Versagene RNA purification kit can isolate total RNA from 3 mL human whole blood collected in Tempus blood RNA tube. Expected yield is 15 to 25 µg RNA per Tempus tube.

Material

1. Tempus blood RNA tube (ABI, 4342792, ABI)
2. Versagene RNA purification kit (Gentra, VGR-0050BP)
(5 PRIME perfectPure RNA Blood Purification Kit, Fisher Scientific)
3. 96-100% ethanol

Method

After blood collection, the Tempus blood RNA tubes are mixed vigorously for 10 sec and are stored at -20°C. If tubes are to be kept at temperatures below -20°C, they are first frozen at -20°C for 24 hr and then are transferred to -80°C.

- equilibrate at room temperature for at least 2 hr (22°C ↓)
 - perform all steps, including centrifugation, at room temperature unless specified
 - prepare DNase Solution
add 2.6 mL DNase Buffer to the lyophilized DNase Enzyme, invert gently to mix, aliquot, and store at -20°C or -80°C
1. after mixing, pour blood from Tempus tube into a clean 50 mL conical tube
 2. add 3 mL 95% ethanol and vortex on high for 2 min
 3. centrifuge at 3000 g for 30 min at 0°C
 4. discard the supernatant
blot on clean absorbent paper and invert tube on absorbent paper for 2 min
 5. add 300 µL Lysis solution and vortex on high speed for 60 sec to thoroughly dissolve the pellet
 6. pipet the lysate onto a purification column
centrifuge at 13000-16000 g for 1 min
transfer the purification column to a new collection tube
 7. add 400 µL Wash 1 Solution and wash the side of the basket during the addition of this solution

- centrifuge at 13000-16000 g for 2 min
transfer carefully the purification column to a new collection tube
8. add 50 μ L DNase Solution and incubate at room temperature for 15 min
 9. add 200 μ L DNase Wash Solution
centrifuge at 13000-16000 g for 1 min
 10. add 200 μ L DNase Wash Solution
centrifuge at 13000-16000 g for 1 min
transfer the purification column to a new tube
 11. add 200 μ L Wash 2 Solution
centrifuge at 13000-16000 g for 1 min
 12. add 200 μ L Wash 2 Solution
centrifuge at 13000-16000 g for 2 min
transfer carefully the purification column to a new tube
if Wash 2 Solution comes in contact with column, centrifuge at
13000-16000 g for 2 min
 13. add 50 μ L Elution Solution and incubate at room temperature for 3 min
centrifuge at 13000-16000 g for 1 min
 14. place the purified RNA on ice and then store at -80°C